



# NovoSci™

## Pediatric Tubing Organizer

NovoSci is pleased to announce the introduction of the Pediatric Tubing Organizer.

The Pediatric Tubing Organizer is ideally suited for

- Pediatric and infant cardiac surgery
- ECMO and post-operative cardiac support
- Anesthesia and Intensive care medicine
- Transplantation teams

### FEATURES

- Easy attachment via secure clamps or suture holes
- Small, compact size with excellent visualization of all lines
- Curved opening to accommodate cautery lines and suction tubing
- Tubing is locked into the holder by the shape of the opening

The NovoSci Pediatric Tubing Organizer will accept the following tubing sizes:

(2) 3/8" lines • (4) 1/4" lines • (2) 3/32" IV lines • (2) 1/8" lines

**Model: 311034-000**

Packaged Single Sterile; 10 per box

**Model: 311034NS**

Packaged Bulk Non-Sterile; 100 per box

Material: Polycarbonate, Stainless Steel • Size: 3.5" x 4"

innovators  
of tomorrow's  
technology  
today



# NovoSci™

2828 North Crescent Ridge Drive  
The Woodlands, TX 77381  
Phone: 281.363.4949  
info@NovoSci.com  
Customer Service: 800.854.0567  
Customer Service Fax: 888.570.4009

Please visit us on the web at <http://www.NovoSci.com>

NovoSci is an FDA registered manufacturer and is ISO 13485:2003 certified